

Product	GUTTA-CUTTER - EASYCUT	
Reference	PRODUCT	
	REFERENCES	
	GUTTA CUTTER EASYCUT	3107172
	GUTTA CUTTER EASYCUT TIPS (4)	3107173
	GUTTA-CUTTER BASE 3TECH	3112377
Picture		
Description	This instrument is meant for the sectioning of gutta-percha points.	
Characteristics & Technical data	<ul style="list-style-type: none"> • Handpiece (rechargeable Ni-MH battery) • Autonomy before recharging: 80 uses • Dimensions of the handpiece: length: 12.6 cm (without cutting tip), \varnothing: 1.9 cm • Mains voltage: 230 V / 50 Hz 	
Storage	Box containing: <ul style="list-style-type: none"> • Instructions for use • A Gutta-Cutter handpiece • A Gutta-Cutter charging stand • 4 interchangeable tips • A Power Supply 	
Indications of use	Cutting tips can be autoclaved up to 25 times. The duration of the sterilization: <ul style="list-style-type: none"> • at 121 ° C at least 20 minutes • at 132/134 ° C at least 5 minutes Always use a rubber dental dam. This instrument is not suitable for use in an environment containing flammable gases.	
Shelf life of the product	12 months.	
Legal Information	Class I medical device. Manufactured by CFPM.	