













MGF HYBRID A 1

Cod. AS-094-EX-AS

Semi-dry 1.2 Hp suction pump, maximum flow rate 1,100 l / min, unerpessure 180 mbar. Integrated motor of secretion separator and amalgam recovery unit designed for 2 workstations. The high reliability of the centrifugal separation system guarantees the correct operation of the suction pump for several years. The pump impellers are made of aluminum and guarantee maximum reliability over time.

-  Maximum simultaneous use: 2
-  Dental Chairs: 2

								
Air capacity		Motor	Voltage Hz	Noise Level	Duty Cycle	max Oper. Head	Size	net Weight
m ³ /min	l/min	kW	230V 50Hz	db(A)	100 %	mbar	mm	Kg
0.0011	1.1	0.94		57		180	570x770x400	22

Maintenance kits

Kit.	Description	Type	Hours
AS-50060240	<p>HYBRID 1, HYBRID 2 annual maintenance kit</p> <p>Comprehends:</p> <ul style="list-style-type: none"> - Water collector - AS-50060214, 1 quantity - Non-return valve, water separator - AS-50060220, 1 quantity - 1.2 m of PVC pipe for water separation - Information Sheet 	Ordinary	1y
AS-50060210	<p>Hybrid 1,2,5 extraordinary maintenance kit</p> <p>The extraordinary maintenance kit includes: liquid separator impeller, non-return valve, 15 mm fitting (with O-ring), clip and information sheet.</p>	Extraordinary	4.000h/4y