

Designation	Powder-free latex finger cots																																		
Reference(s)	Size	S	M	L	XL																														
	Reference	3160736	3160737	3160738	3160739																														
	UDI-DI	06958627514209	06958627514216	06958627514223	06958627514230																														
Visual																																			
Description	Single-use, non-sterile, powder-free finger cot to protect the finger, probe or sensor during medical examinations.																																		
Specifications and technical data	<table><tr><th>Sizes</th><th>S</th><th>M</th><th>L</th><th>XL</th><th>Tolerance</th></tr><tr><td>Width (mm)</td><td>24</td><td>28</td><td>30</td><td>34.5</td><td>+/- 2 mm</td></tr><tr><td>Length (mm)</td><td colspan="4">70</td><td>+/- 5 mm</td></tr><tr><td>Thickness (mm)</td><td colspan="4">0.09</td><td>+/- 0.03 mm</td></tr><tr><td>Average weight (g)</td><td>0.41</td><td>0.46</td><td>0.52</td><td>0.58</td><td>+/- 0.02 g</td></tr></table>					Sizes	S	M	L	XL	Tolerance	Width (mm)	24	28	30	34.5	+/- 2 mm	Length (mm)	70				+/- 5 mm	Thickness (mm)	0.09				+/- 0.03 mm	Average weight (g)	0.41	0.46	0.52	0.58	+/- 0.02 g
	Sizes	S	M	L	XL	Tolerance																													
	Width (mm)	24	28	30	34.5	+/- 2 mm																													
	Length (mm)	70				+/- 5 mm																													
	Thickness (mm)	0.09				+/- 0.03 mm																													
Average weight (g)	0.41	0.46	0.52	0.58	+/- 0.02 g																														
	<table><tr><th>Test</th><th>Spéc.</th><th>Reference standard</th></tr><tr><td>Tensile strength (Mpa)</td><td>≥ 24</td><td rowspan="2">ASTM D3772/D412</td></tr><tr><td>Elongation at break (%)</td><td>≥ 750</td></tr></table>					Test	Spéc.	Reference standard	Tensile strength (Mpa)	≥ 24	ASTM D3772/D412	Elongation at break (%)	≥ 750																						
Test	Spéc.	Reference standard																																	
Tensile strength (Mpa)	≥ 24	ASTM D3772/D412																																	
Elongation at break (%)	≥ 750																																		
Packaging	Bag of 100 finger cots carton of 100 bags																																		
Shelf life	5 years																																		
Indications for use																																			
Legal information	Class I medical device under MDR regulation 2017/745 CE mark - Manufacturer: Ningbo Greetmed Medical Instruments Co, Ltd.																																		

Version 1.0

Date : Updated 31/10/2023 by LF